

ERCC6L2 Polyclonal Antibody

Catalog No: #27919

Package Size: #27919-1 50ul #27919-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	ERCC6L2 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Hu,Ms
Immunogen Description	Recombinant fusion protein of human ERCC6L2 (NP_064592.2).
Other Names	ERCC6L2; BMFS2; C9orf102; RAD26L; SR278; ERCC excision repair 6 like 2
Accession No.	Swiss-Prot#:Q5T890NCBI Gene ID:375748
Calculated MW	81kDa
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

Application Details

WB □ 1:500 - 1:2000

Images



Western blot analysis of extracts of various cell lines, using ERCC6L2 antibody.

Background

This gene encodes a member of the Snf2 family of helicase-like proteins. The encoded protein may play a role in DNA repair and mitochondrial function. Mutations in this gene have been associated with bone marrow failure syndrome 2. Alternatively spliced transcript variants that encode different protein isoforms have been described.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.