

# PID1 Polyclonal Antibody

Catalog No: #27951

Package Size: #27951-1 50ul #27951-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)


## Description

Product Name	PID1 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Hu,Ms,Rt
Immunogen Description	Recombinant fusion protein of human PID1 (NP_060403.3).
Other Names	PID1; HMFN2073; NYGGF4; P-CL11; PCL11; phosphotyrosine interaction domain containing 1
Accession No.	Swiss-Prot#:Q7Z2X4NCBI Gene ID:55022
Calculated MW	27kDa
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

## Application Details

WB □ 1:500 - 1:2000

## Images



Western blot analysis of extracts of various cell lines, using PID1 antibody.

## Background

## Published Papers

el et al., Dual role of PID1 in regulating apoptosis induced by distinct anticancer-agents through AKT/Raf-1-dependent pathway in hepatocellular carcinomaInCell Death DiscovOn2023 Apr 28byJian Yang,Senlin Li et al..PMID:37117198, , (2023)

[PMID:37117198](https://pubmed.ncbi.nlm.nih.gov/37117198/)

---

Note: This product is for in vitro research use only and is not intended for use in humans or animals.