

DBNL Polyclonal Antibody

Catalog No: #28333



Package Size: #28333-1 50ul #28333-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	DBNL Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB,IF
Species Reactivity	Hu,Ms,Rt
Immunogen Description	Recombinant fusion protein of human DBNL (NP_001014436.1).
Other Names	DBNL; ABP1; HIP-55; HIP55; SH3P7; drebrin like
Accession No.	Swiss-Prot#:Q9UJU6NCBI Gene ID:28988
Calculated MW	56kDa
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

Application Details

WB 1:500 - 1:2000 IF 1:50 - 1:100

Images



Western blot analysis of extracts of various cell lines, using DBNL antibody.

Background

Note: This product is for in vitro research use only and is not intended for use in humans or animals.