

# ZNF396 Polyclonal Antibody

Catalog No: #28372



Package Size: #28372-1 50ul #28372-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)


## Description

Product Name	ZNF396 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Hu,Ms,Rt
Immunogen Description	Recombinant fusion protein of human ZNF396 (NP_665699.1).
Other Names	ZNF396; ZSCAN14; zinc finger protein 396
Accession No.	Swiss-Prot#:Q96N95NCBI Gene ID:252884
Calculated MW	41kDa
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

## Application Details

WB □ 1:500 - 1:2000

## Images



Western blot analysis of extracts of various cell lines, using ZNF396 antibody.

## Background

Note: This product is for in vitro research use only and is not intended for use in humans or animals.