

## ZNF333 Polyclonal Antibody

Catalog No: #28381

Package Size: #28381-1 50ul #28381-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


## Description

|                       |   |
|-----------------------|---|
| Product Name          | ZNF333 Polyclonal Antibody  |
| Host Species          | Rabbit  |
| Clonality             | Polyclonal  |
| Isotype               | IgG   |
| Purification          | Affinity purification   |
| Applications          | WB  |
| Species Reactivity    | Hu,Ms,Rt  |
| Immunogen Description | Recombinant fusion protein of human ZNF333 (NP_001287841.1).        |
| Other Names           | ZNF333; zinc finger protein 333                                     |
| Accession No.         | Swiss-Prot#:Q96JL9NCBI Gene ID:84449                                |
| Calculated MW         | 80kDa   |
| Formulation           | Avoid freeze / thaw cycles.   Buffer: PBS with 50% glycerol, pH7.4. |
| Storage               | Store at -20°C  |

## Application Details

WB □ 1:500 - 1:2000

## Images



Western blot analysis of extracts of various cell lines, using ZNF333 antibody.

## Background

Note: This product is for in vitro research use only and is not intended for use in humans or animals.