

## DMKN Polyclonal Antibody

Catalog No: #28391

Package Size: #28391-1 50ul #28391-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)


## Description

Product Name	DMKN Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Hu,Rt
Immunogen Description	Recombinant fusion protein of human DMKN (NP_001177276.1).
Other Names	DMKN; UNQ729; ZD52F10; dermokine
Accession No.	Swiss-Prot#:Q6E0U4NCBI Gene ID:93099
Calculated MW	47kDa
Formulation	Avoid freeze / thaw cycles.   Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

## Application Details

WB □ 1:500 - 1:2000

## Images



Western blot analysis of extracts of rat ovary, using DMKN antibody.

## Background

This gene is upregulated in inflammatory diseases, and it was first observed as expressed in the differentiated layers of skin. The most interesting aspect of this gene is the differential use of promoters and terminators to generate isoforms with unique cellular distributions and domain components. Alternatively spliced transcript variants encoding different isoforms have been identified for this gene.

---

Note: This product is for in vitro research use only and is not intended for use in humans or animals.