

ELMOD2 Polyclonal Antibody

Catalog No: #28457

Package Size: #28457-1 50ul #28457-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	ELMOD2 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Mouse
Immunogen Description	Recombinant fusion protein of human ELMOD2 (NP_714913.1).
Other Names	ELMOD2; 9830169G11Rik; ELMO domain containing 2
Accession No.	Swiss-Prot#:Q8IZ81NCBI Gene ID:255520
Calculated MW	35kDa
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

Application Details

WB □ 1:500 - 1:2000

Images



Western blot analysis of extracts of various cell lines, using ELMOD2 antibody at 1:500 dilution.

Background

This gene encodes one of six engulfment and motility (ELMO) domain-containing proteins. This gene is thought to play a role in antiviral responses. Mutations in this gene may be involved in the cause of familial idiopathic pulmonary fibrosis.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.