

ZFR Polyclonal Antibody

Catalog No: #28479

Package Size: #28479-1 50ul #28479-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

| | |
|-----------------------|---|
| Product Name | ZFR Polyclonal Antibody |
| Host Species | Rabbit |
| Clonality | Polyclonal |
| Isotype | IgG |
| Purification | Affinity purification |
| Applications | WB |
| Species Reactivity | Hu,Ms,Rt |
| Immunogen Description | Recombinant fusion protein of human ZFR (NP_057191.2). |
| Other Names | ZFR; SPG71; ZFR1; zinc finger RNA-binding protein |
| Accession No. | Swiss-Prot#:Q96KR1NCBI Gene ID:51663 |
| Calculated MW | 150kDa |
| Formulation | Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4. |
| Storage | Store at -20°C |

Application Details

WB □ 1:500 - 1:2000

Images



Western blot analysis of extracts of various cell lines, using ZFR antibody.

Background

Note: This product is for in vitro research use only and is not intended for use in humans or animals.