

ZSCAN31 Polyclonal Antibody

Catalog No: #28486

Package Size: #28486-1 50ul #28486-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

| | |
|-----------------------|---|
| Product Name | ZSCAN31 Polyclonal Antibody |
| Host Species | Rabbit |
| Clonality | Polyclonal |
| Isotype | IgG |
| Purification | Affinity purification |
| Applications | WB |
| Species Reactivity | Rt |
| Immunogen Description | Recombinant fusion protein of human ZSCAN31 (NP_112161.3). |
| Other Names | ZSCAN31; ZNF20-Lp; ZNF310P; ZNF323; zinc finger and SCAN domain containing 31 |
| Accession No. | Swiss-Prot#:Q96LW9NCBI Gene ID:64288 |
| Calculated MW | 47kDa |
| Formulation | Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4. |
| Storage | Store at -20°C |

Application Details

WB □ 1:500 - 1:2000

Images



Western blot analysis of extracts of rat liver, using ZSCAN31 antibody.

Background

This gene encodes a protein containing multiple C2H2-type zinc finger motifs. Alternative splicing results in multiple transcript variants.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.