

ZNF319 Polyclonal Antibody

Catalog No: #28506

Package Size: #28506-1 50ul #28506-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	ZNF319 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Hu,Ms,Rt
Immunogen Description	Recombinant fusion protein of human ZNF319 (NP_065858.1).
Other Names	ZNF319; ZFP319; zinc finger protein 319
Accession No.	Swiss-Prot#:Q9P2F9NCBI Gene ID:57567
Calculated MW	65kDa
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

Application Details

WB □ 1:500 - 1:2000

Images



Western blot analysis of extracts of various cell lines, using ZNF319 antibody.

Background

Note: This product is for in vitro research use only and is not intended for use in humans or animals.