

ZNF101 Polyclonal Antibody

Catalog No: #28526

Package Size: #28526-1 50ul #28526-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	ZNF101 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Human,Mouse,Rat
Immunogen Description	Recombinant fusion protein of human ZNF101 (NP_149981.2).
Other Names	ZNF101; HZF12; zinc finger protein 101
Accession No.	Swiss-Prot#:Q8IZC7NCBI Gene ID:94039
Calculated MW	50kDa
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

Application Details

WB □ 1:500 - 1:2000

Images



Western blot analysis of extracts of various cell lines, using ZNF101 antibody at 1:1000 dilution.

Background

Note: This product is for in vitro research use only and is not intended for use in humans or animals.