

ZSCAN20 Polyclonal Antibody

Catalog No: #28543

Package Size: #28543-1 50ul #28543-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	ZSCAN20 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Hu,Ms
Immunogen Description	Recombinant fusion protein of human ZSCAN20 (NP_660281.2).
Other Names	ZSCAN20; KOX29; ZFP-31; ZNF31; ZNF360; zinc finger and SCAN domain containing 20
Accession No.	Swiss-Prot#:P17040NCBI Gene ID:7579
Calculated MW	117kDa
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

Application Details

WB □ 1:500 - 1:2000

Images



Western blot analysis of extracts of various cell lines, using ZSCAN20 antibody.

Background

Note: This product is for in vitro research use only and is not intended for use in humans or animals.