

ZNF699 Polyclonal Antibody

Catalog No: #28598

Package Size: #28598-1 50ul #28598-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	ZNF699 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Hu
Immunogen Description	Recombinant fusion protein of human ZNF699 (NP_940937.1).
Other Names	ZNF699; hang; zinc finger protein 699
Accession No.	Swiss-Prot#:Q32M78NCBI Gene ID:374879
Calculated MW	74kDa
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

Application Details

WB □ 1:500 - 1:2000

Images



Western blot analysis of extracts of various cell lines, using ZNF699 antibody.

Background

Note: This product is for in vitro research use only and is not intended for use in humans or animals.