

DCUN1D1 Polyclonal Antibody

Catalog No: #28633

Package Size: #28633-1 50ul #28633-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	DCUN1D1 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Hu,Ms
Immunogen Description	Recombinant fusion protein of human DCUN1D1 (NP_065691.2).
Other Names	DCUN1D1; DCNL1; DCUN1L1; RP42; SCCRO; SCRO; Tes3; DCN1-like protein 1
Accession No.	Swiss-Prot#:Q96GG9NCBI Gene ID:54165
Calculated MW	30kDa
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

Application Details

WB □ 1:500 - 1:2000

Images



Western blot analysis of extracts of various cell lines, using DCUN1D1 antibody.

Background

Note: This product is for in vitro research use only and is not intended for use in humans or animals.