

DCSTAMP Polyclonal Antibody

Catalog No: #28655

Package Size: #28655-1 50ul #28655-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	DCSTAMP Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Human,Mouse
Immunogen Description	Recombinant fusion protein of human DCSTAMP (NP_110415.1).
Other Names	DCSTAMP;FIND;TM7SF4;hDC-STAMP
Accession No.	Uniprot:Q9H295GeneID:81501
Calculated MW	53kDa
SDS-PAGE MW	53KDa
Formulation	PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.

Application Details

WB □ 1:500 - 1:2000

Images



Western blot analysis of extracts of Mouse brain, using DCSTAMP antibody.

Background

This gene encodes a seven-pass transmembrane protein that is primarily expressed in dendritic cells. The encoded protein is involved in a range of immunological functions carried out by dendritic cells. This protein plays a role in osteoclastogenesis and myeloid differentiation. Alternate splicing results in multiple transcript variants.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.