

## ZNF296 Polyclonal Antibody

Catalog No: #28657

Package Size: #28657-1 50ul #28657-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


## Description

Product Name	ZNF296 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Hu,Ms
Immunogen Description	Recombinant fusion protein of human ZNF296 (NP_660331.1).
Other Names	ZNF296; ZFP296; ZNF342; zinc finger protein 296
Accession No.	Swiss-Prot#:Q8WUU4NCBI Gene ID:162979
Calculated MW	56kDa
Formulation	Avoid freeze / thaw cycles.   Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

## Application Details

WB □ 1:500 - 1:2000

## Images



Western blot analysis of extracts of various cell lines, using ZNF296 antibody.

## Background

Note: This product is for in vitro research use only and is not intended for use in humans or animals.