

ZNF141 Polyclonal Antibody

Catalog No: #28707

Package Size: #28707-1 50ul #28707-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	ZNF141 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	IHC,IF
Species Reactivity	Hu,Ms,Rt
Immunogen Description	Recombinant fusion protein of human ZNF141 (NP_003432.1).
Other Names	ZNF141; D4S90; pHZ-44; PAPA6; zinc finger protein 141
Accession No.	Swiss-Prot#:Q15928NCBI Gene ID:7700
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

Application Details

IHC□1:50 - 1:200IF□1:50 - 1:200

Images



Immunohistochemistry of paraffin-embedded mouse lung using ZNF141 antibody.

Background

Note: This product is for in vitro research use only and is not intended for use in humans or animals.