

ZNF177 Polyclonal Antibody

Catalog No: #28709

Package Size: #28709-1 50ul #28709-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

| | |
|-----------------------|---|
| Product Name | ZNF177 Polyclonal Antibody |
| Host Species | Rabbit |
| Clonality | Polyclonal |
| Isotype | IgG |
| Purification | Affinity purification |
| Applications | WB |
| Species Reactivity | Ms |
| Immunogen Description | Recombinant fusion protein of human ZNF177 (NP_003442.2). |
| Other Names | ZNF177; PIGX; zinc finger protein 177 |
| Accession No. | Swiss-Prot#:Q13360NCBI Gene ID:7730 |
| Calculated MW | 55kDa |
| Formulation | Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4. |
| Storage | Store at -20°C |

Application Details

WB □ 1:500 - 1:2000

Images



Western blot analysis of extracts of mouse testis, using ZNF177 antibody.

Background

Note: This product is for in vitro research use only and is not intended for use in humans or animals.