

ELP5 Polyclonal Antibody

Catalog No: #28732

Package Size: #28732-1 50ul #28732-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	ELP5 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Hu,Ms
Immunogen Description	Recombinant fusion protein of human ELP5 (NP_056177.3).
Other Names	ELP5; C17orf81; DERP6; HSPC002; MST071; MSTP071; elongator complex protein 5
Accession No.	Swiss-Prot#:Q8TE02NCBI Gene ID:23587
Calculated MW	38kDa
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

Application Details

WB □ 1:500 - 1:2000

Images



Western blot analysis of extracts of various cell lines, using ELP5 antibody.

Background

Note: This product is for in vitro research use only and is not intended for use in humans or animals.