

DCAF4 Polyclonal Antibody

Catalog No: #28735

Package Size: #28735-1 50ul #28735-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	DCAF4 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Hu
Immunogen Description	Recombinant fusion protein of human DCAF4 (NP_056419.2).
Other Names	DCAF4; WDR21; WDR21A; DDB1 and CUL4 associated factor 4
Accession No.	Swiss-Prot#:Q8WV16NCBI Gene ID:26094
Calculated MW	56kDa
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

Application Details

WB □ 1:500 - 1:2000

Images



Western blot analysis of extracts of U-251MG cells, using DCAF4 antibody.

Background

This gene encodes a WD repeat-containing protein that interacts with the Cul4-Ddb1 E3 ligase macromolecular complex. Multiple alternatively spliced transcript variants have been found for this gene.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.