

ZNF415 Polyclonal Antibody

Catalog No: #28755



Package Size: #28755-1 50ul #28755-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	ZNF415 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB,IHC
Species Reactivity	Hu,Ms,Rt
Immunogen Description	Recombinant fusion protein of human ZNF415 (NP_060825.2).
Other Names	ZNF415; Zfp; zinc finger protein 415
Accession No.	Swiss-Prot#:Q09FC8NCBI Gene ID:55786
Calculated MW	Refer to figures
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C


Application Details

WB 1:500 - 1:2000 IHC 1:50 - 1:200


Images



Immunohistochemistry of paraffin-embedded rat brain using ZNF415 antibody.



Immunohistochemistry of paraffin-embedded mouse brain using ZNF415 antibody.



Immunohistochemistry of paraffin-embedded mouse kidney using ZNF415 antibody.

Background

Note: This product is for in vitro research use only and is not intended for use in humans or animals.