

VIPAS39 Polyclonal Antibody

Catalog No: #28757

Package Size: #28757-1 50ul #28757-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	VIPAS39 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB,IHC
Species Reactivity	Ms,Rt
Immunogen Description	Recombinant fusion protein of human VIPAS39 (NP_071350.2).
Other Names	VIPAS39; C14orf133; SPE-39; SPE39; VIPAR; VPS16B; hSPE-39; spermatogenesis-defective protein 39 homolog
Accession No.	Swiss-Prot#:Q9H9C1NCBI Gene ID:63894
Calculated MW	Refer to figures
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

Application Details

WB □ 1:500 - 1:2000 IHC □ 1:50 - 1:200

Background

This gene encodes a protein involved in the sorting of lysosomal proteins. Mutations in this gene are associated with ARCS2 (arthrogryposis, renal dysfunction, and cholestasis-2). Alternative splicing results in multiple transcript variants.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.