

ZNF549 Polyclonal Antibody

Catalog No: #28790

Package Size: #28790-1 50ul #28790-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

| | |
|-----------------------|---|
| Product Name | ZNF549 Polyclonal Antibody |
| Host Species | Rabbit |
| Clonality | Polyclonal |
| Isotype | IgG |
| Purification | Affinity purification |
| Applications | WB,IHC |
| Species Reactivity | Hu |
| Immunogen Description | Recombinant fusion protein of human ZNF549 (NP_001186224.1). |
| Other Names | ZNF549; zinc finger protein 549 |
| Accession No. | Swiss-Prot#:Q6P9A3NCBI Gene ID:256051 |
| Calculated MW | 68kDa |
| Formulation | Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4. |
| Storage | Store at -20°C |

Application Details

WB □ 1:500 - 1:2000 IHC □ 1:50 - 1:200

Images



Western blot analysis of extracts of Human placenta, using ZNF549 antibody.

Background

Note: This product is for in vitro research use only and is not intended for use in humans or animals.