

ZNF677 Polyclonal Antibody

Catalog No: #28792

Package Size: #28792-1 50ul #28792-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

| | |
|-----------------------|---|
| Product Name | ZNF677 Polyclonal Antibody |
| Host Species | Rabbit |
| Clonality | Polyclonal |
| Isotype | IgG |
| Purification | Affinity purification |
| Applications | WB |
| Species Reactivity | Rt |
| Immunogen Description | Recombinant fusion protein of human ZNF677 (NP_872415.1). |
| Other Names | ZNF677; zinc finger protein 677 |
| Accession No. | Swiss-Prot#:Q86XU0NCBI Gene ID:342926 |
| Calculated MW | 68kDa |
| Formulation | Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4. |
| Storage | Store at -20°C |

Application Details

WB □ 1:500 - 1:2000

Images



Western blot analysis of extracts of rat testis, using ZNF677 antibody.

Background

Note: This product is for in vitro research use only and is not intended for use in humans or animals.