

## QSOX2 Polyclonal Antibody

Catalog No: #28873

Package Size: #28873-1 50ul #28873-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)


## Description

Product Name	QSOX2 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Ms
Immunogen Description	Recombinant fusion protein of human QSOX2 (NP_859052.3).
Other Names	QSOX2; QSCN6L1; SOXN; sulfhydryl oxidase 2
Accession No.	Swiss-Prot#:Q6ZRP7NCBI Gene ID:169714
Calculated MW	85kDa
Formulation	Avoid freeze / thaw cycles.   Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

## Application Details

WB □ 1:500 - 1:2000

## Images



Western blot analysis of extracts of various cell lines, using QSOX2 antibody.

## Background

Note: This product is for in vitro research use only and is not intended for use in humans or animals.