

## RBPMS Polyclonal Antibody

Catalog No: #28946

Package Size: #28946-1 50ul #28946-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)


## Description

Product Name	RBPMS Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Hu,Ms,Rt
Immunogen Description	Recombinant fusion protein of human RBPMS (NP_001008712.1).
Other Names	RBPMS; HERMES; RNA-binding protein with multiple splicing
Accession No.	Swiss-Prot#:Q93062NCBI Gene ID:11030
Calculated MW	22kDa
Formulation	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.

## Application Details

WB □ 1:200 - 1:2000

## Images



Western blot analysis of RBPMS expression in mouse liver (A), mouse lung (B), rat testis (C) whole cell lysates.

## Background

This gene encodes a member of the RNA recognition motif family of RNA-binding proteins. The RNA recognition motif is between 80-100 amino acids in length and family members contain one to four copies of the motif. The RNA recognition motif consists of two short stretches of conserved sequence, as well as a few highly conserved hydrophobic residues. The encoded protein has a single, putative RNA recognition motif in its N-terminus. Alternative splicing results in multiple transcript variants encoding different isoforms.

---

Note: This product is for in vitro research use only and is not intended for use in humans or animals.