

ESYT1 Polyclonal Antibody

Catalog No: #28954



Package Size: #28954-1 50ul #28954-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	ESYT1 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB,IF
Species Reactivity	Human,Mouse,Rat
Immunogen Description	Recombinant fusion protein of human ESYT1 (NP_056107.1).
Other Names	ESYT1; FAM62A; MBC2; extended synaptotagmin-1
Accession No.	Swiss-Prot#:Q9BSJ8NCBI Gene ID:23344
Calculated MW	123kDa
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C


Application Details

WB □ 1:200 - 1:2000 IF □ 1:50 - 1:200

Images



Immunofluorescence analysis of U-2 OS cells using ESYT1 antibody at dilution of 1:100. Blue: DAPI for nuclear staining.



Western blot analysis of extracts of various cell lines, using ESYT1 antibody at 1:1000 dilution.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.