

DRAM1 Polyclonal Antibody

Catalog No: #28982

Package Size: #28982-1 50ul #28982-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	DRAM1 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Ms
Immunogen Description	A synthetic peptide of human DRAM1 (NP_060840.2).
Other Names	DRAM1; DRAM; DNA damage regulated autophagy modulator 1
Accession No.	Swiss-Prot#:Q8N682NCBI Gene ID:55332
Calculated MW	26kDa
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

Application Details

WB □ 1:200 - 1:2000

Images



Western blot analysis of extracts of various cell lines, using DRAM1 antibody.

Background

This gene is regulated as part of the p53 tumor suppressor pathway. The gene encodes a lysosomal membrane protein that is required for the induction of autophagy by the pathway. Decreased transcriptional expression of this gene is associated with various tumors. This gene has a pseudogene on chromosome 4.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.