

## ZNF331 Polyclonal Antibody

Catalog No: #28984

Package Size: #28984-1 50ul #28984-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)


## Description

Product Name	ZNF331 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Ms,Rt
Immunogen Description	Recombinant fusion protein of human ZNF331 (NP_061025.5).
Other Names	ZNF331; RITA; ZNF361; ZNF463; zinc finger protein 331
Accession No.	Swiss-Prot#:Q9NQX6NCBI Gene ID:55422
Calculated MW	56kDa
Formulation	Avoid freeze / thaw cycles.   Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

## Application Details

WB □ 1:200 - 1:2000

## Images



Western blot analysis of extracts of various cell lines, using ZNF331 antibody.

## Background

This gene encodes a zinc finger protein containing a KRAB (Kruppel-associated box) domain found in transcriptional repressors. This gene may be methylated and silenced in cancer cells. This gene is located within a differentially methylated region (DMR) and shows allele-specific expression in placenta. Alternative splicing and the use of alternative promoters results in multiple transcript variants encoding the same protein.

---

Note: This product is for in vitro research use only and is not intended for use in humans or animals.