

Zfp541 Polyclonal Antibody

Catalog No: #29324

Package Size: #29324-1 50ul #29324-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	Zfp541 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Ms,Rt
Immunogen Description	Recombinant fusion protein of mouse Zfp541 (NP_001092747.1).
Other Names	ZNF541; zinc finger protein 541
Accession No.	Swiss-Prot#:_NCBI Gene ID:666528
Calculated MW	170kDa
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

Application Details

WB 1:200 - 1:2000

Images



Western blot analysis of extracts of various cell lines, using ZNF541 antibody.

Background

Note: This product is for in vitro research use only and is not intended for use in humans or animals.