

ZNF502 Polyclonal Antibody

Catalog No: #29334

Package Size: #29334-1 50ul #29334-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	ZNF502 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Hu
Immunogen Description	Recombinant fusion protein of human ZNF502 (NP_149987.2).
Other Names	ZNF502; zinc finger protein 502
Accession No.	Swiss-Prot#:Q8TBZ5NCBI Gene ID:91392
Calculated MW	70kDa
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

Application Details

WB □ 1:200 - 1:2000

Images



Western blot analysis of extracts of 293T cells, using ZNF502 antibody.

Background

Note: This product is for in vitro research use only and is not intended for use in humans or animals.