

DNTTIP1 Polyclonal Antibody

Catalog No: #29339



Package Size: #29339-1 50ul #29339-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	DNTTIP1 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Hu,Ms,Rt
Immunogen Description	Recombinant fusion protein of human DNTTIP1 (NP_443183.1).
Other Names	DNTTIP1; C20orf167; Tdif1; dJ447F3.4; deoxynucleotidyltransferase terminal-interacting protein 1
Accession No.	Swiss-Prot#:Q9H147NCBI Gene ID:116092
Calculated MW	37kDa
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

Application Details

WB □ 1:200 - 1:2000

Images



Western blot analysis of extracts of various cell lines, using DNTTIP1 antibody.

Background

Note: This product is for in vitro research use only and is not intended for use in humans or animals.