

DNAJC24 Polyclonal Antibody

Catalog No: #29340

Package Size: #29340-1 50ul #29340-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	DNAJC24 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Ms
Immunogen Description	Recombinant fusion protein of human DNAJC24 (NP_859057.4).
Other Names	DNAJC24; DPH4; JJJ3; ZCSL3; dnaJ homolog subfamily C member 24
Accession No.	Swiss-Prot#:Q6P3W2NCBI Gene ID:120526
Calculated MW	17kDa
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

Application Details

WB □ 1:200 - 1:2000

Images



Western blot analysis of extracts of various cell lines, using DNAJC24 antibody.

Background

Note: This product is for in vitro research use only and is not intended for use in humans or animals.