

ZNF133 Polyclonal Antibody

Catalog No: #29404

Package Size: #29404-1 50ul #29404-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	ZNF133 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Hu,Ms,Rt
Immunogen Description	Recombinant fusion protein of human ZNF133 (NP_003425.2).
Other Names	ZNF133; ZNF150; pHZ-13; pHZ-66; zinc finger protein 133
Accession No.	Swiss-Prot#:P52736NCBI Gene ID:7692
Calculated MW	70kDa
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

Application Details

WB □ 1:500 - 1:2000

Images



Western blot analysis of extracts of various cell lines, using ZNF133 antibody.

Background

Note: This product is for in vitro research use only and is not intended for use in humans or animals.