

Defb26 Polyclonal Antibody

Catalog No: #29483

Package Size: #29483-1 50ul #29483-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

| | |
|-----------------------|---|
| Product Name | Defb26 Polyclonal Antibody |
| Host Species | Rabbit |
| Clonality | Polyclonal |
| Isotype | IgG |
| Purification | Affinity purification |
| Applications | IF |
| Species Reactivity | Hu |
| Immunogen Description | Recombinant fusion protein of rat Defb26 (NP_001032595.1). |
| Other Names | DEFB126; C20orf8; DEFB-26; DEFB26; HBD26; bA530N10.1; hBD-26; beta-defensin 126 |
| Accession No. | Swiss-Prot#:Q32ZG6NCBI Gene ID:641629 |
| Formulation | Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4. |
| Storage | Store at -20°C |

Application Details

IF □ 1:50 - 1:200

Images



Immunofluorescence analysis of U2OS cells using Defb26 antibody.

Background

Note: This product is for in vitro research use only and is not intended for use in humans or animals.