

## ZNF517 Polyclonal Antibody

Catalog No: #29524

Package Size: #29524-1 50ul #29524-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)


## Description

Product Name	ZNF517 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Hu
Immunogen Description	Recombinant fusion protein of human ZNF517 (NP_998770.2).
Other Names	ZNF517; zinc finger protein 517
Accession No.	Swiss-Prot#:Q6ZMY9NCBI Gene ID:340385
Calculated MW	56kDa
Formulation	Avoid freeze / thaw cycles.   Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

## Application Details

WB □ 1:500 - 1:2000

## Images



Western blot analysis of extracts of K-562 cells, using ZNF517 antibody.

## Background

Note: This product is for in vitro research use only and is not intended for use in humans or animals.