

ZNF644 Polyclonal Antibody

Catalog No: #29754

Package Size: #29754-1 50ul #29754-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	ZNF644 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Ms
Immunogen Description	Recombinant fusion protein of human ZNF644 (NP_958357.1).
Other Names	ZNF644; BM-005; MYP21; NatF; ZEP-2; zinc finger protein 644
Accession No.	Swiss-Prot#:Q9H582NCBI Gene ID:84146
Calculated MW	150kDa
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

Application Details

WB □ 1:500 - 1:2000

Images



Western blot analysis of extracts of mouse liver, using ZNF644 antibody.

Background

The protein encoded by this gene is a zinc finger transcription factor that may play a role in eye development. Defects in this gene have been associated with high myopia. Three transcript variants encoding two different isoforms have been found for this gene.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.