

ZNF587 Polyclonal Antibody

Catalog No: #29755

Package Size: #29755-1 50ul #29755-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	ZNF587 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Hu,Ms,Rt
Immunogen Description	Recombinant fusion protein of human ZNF587 (NP_116217.1).
Other Names	ZNF587; ZF6; zinc finger protein 587
Accession No.	Swiss-Prot#:Q96SQ5NCBI Gene ID:84914
Calculated MW	68kDa
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

Application Details

WB □ 1:500 - 1:2000

Images



Western blot analysis of extracts of mouse heart, using ZNF587 antibody.

Background

Note: This product is for in vitro research use only and is not intended for use in humans or animals.