

EFCAB3 Polyclonal Antibody

Catalog No: #29759

Package Size: #29759-1 50ul #29759-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	EFCAB3 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Ms,Rt
Immunogen Description	Recombinant fusion protein of human EFCAB3 (NP_775774.1).
Other Names	EFCAB3; EF-hand calcium binding domain 3
Accession No.	Swiss-Prot#:Q8N7B9NCBI Gene ID:146779
Calculated MW	55kDa
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

Application Details

WB □ 1:500 - 1:2000

Images



Western blot analysis of extracts of various cell lines, using EFCAB3 antibody.

Background

Note: This product is for in vitro research use only and is not intended for use in humans or animals.