

## MAGOHB Polyclonal Antibody

Catalog No: #29770



Package Size: #29770-1 50ul #29770-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)


## Description

Product Name	MAGOHB Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Human,Rat
Immunogen Description	Recombinant fusion protein of human MAGOHB (NP_060518.1).
Other Names	MAGOHB;MGN2;mago;magoh
Accession No.	Uniprot:Q96A72GeneID:55110
Calculated MW	Refer to figures
SDS-PAGE MW	17kDa
Formulation	PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.

## Application Details

WB □ 1:500 - 1:2000

## Images



Western blot analysis of extracts of various cell lines, using MAGOHB Rabbit pAb.

## Background

Required for pre-mRNA splicing as component of the spliceosome. Plays a redundant role with MAGOH in the exon junction complex and in the nonsense-mediated decay (NMD) pathway.

---

Note: This product is for in vitro research use only and is not intended for use in humans or animals.