

## VASN Polyclonal Antibody

Catalog No: #29780

Package Size: #29780-1 50ul #29780-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)


## Description

Product Name	VASN Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Mouse,Rat
Immunogen Description	Recombinant fusion protein of human VASN (NP_612449.2).
Other Names	VASN; SLITL2; vasorin
Accession No.	Swiss-Prot#:Q6EMK4NCBI Gene ID:114990
Calculated MW	71kDa
Formulation	Avoid freeze / thaw cycles.   Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

## Application Details

WB □ 1:500 - 1:2000

## Images



Western blot analysis of extracts of various cell lines, using VASN antibody at 1:1000 dilution.

## Background

Note: This product is for in vitro research use only and is not intended for use in humans or animals.