

## DPPA5 Polyclonal Antibody

Catalog No: #29794

Package Size: #29794-1 50ul #29794-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


## Description

Product Name	DPPA5 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Ms,Rt
Immunogen Description	Recombinant fusion protein of human DPPA5 (NP_001020461.1).
Other Names	DPPA5; ESG1; developmental pluripotency associated 5
Accession No.	Swiss-Prot#:A6NC42NCBI Gene ID:340168
Calculated MW	13kDa
Formulation	Avoid freeze / thaw cycles.   Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

## Application Details

WB □ 1:500 - 1:2000

## Images



Western blot analysis of extracts of various cell lines, using DPPA5 antibody.

## Background

This gene encodes a protein that may function in the control of cell pluripotency and early embryogenesis. Expression of this gene is a specific marker for pluripotent stem cells. Pseudogenes of this gene are located on the short arm of chromosome 10 and the long arm of chromosomes 14 and 19.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.