

## ZNF471 Polyclonal Antibody

Catalog No: #29882

Package Size: #29882-1 50ul #29882-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


## Description

Product Name	ZNF471 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Ms,Rt
Immunogen Description	Recombinant fusion protein of human ZNF471 (NP_065864.2).
Other Names	ZNF471; ERP1; Z1971; zinc finger protein 471
Accession No.	Swiss-Prot#:Q9BX82NCBI Gene ID:57573
Calculated MW	90kDa
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

## Application Details

WB □ 1:500 - 1:2000

## Images



Western blot analysis of extracts of various cell lines, using ZNF471 antibody.

## Background

Note: This product is for in vitro research use only and is not intended for use in humans or animals.