

ARL6IP4 Polyclonal Antibody

Catalog No: #29984

Package Size: #29984-1 50ul #29984-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	ARL6IP4 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	IHC
Species Reactivity	Hu
Immunogen Description	Recombinant fusion protein of human ARL6IP4 (NP_057722.3).
Other Names	ARL6IP4; SFRS20; SR-25; SRp25; SRrp37; ADP ribosylation factor like GTPase 6 interacting protein 4
Accession No.	Swiss-Prot#:Q66PJ3NCBI Gene ID:51329
Calculated MW	Refer to figures
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C


Application Details

IHC □ 1:50 - 1:200

Images



Immunohistochemistry of paraffin-embedded human lung using ARL6IP4 antibody.



Immunohistochemistry of paraffin-embedded human breast cancer using ARL6IP4 antibody.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.