

## ZNF563 Polyclonal Antibody

Catalog No: #30002

Package Size: #30002-1 50ul #30002-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	ZNF563 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Hu,Ms,Rt
Immunogen Description	A synthetic peptide of human ZNF563 (NP_660319.1).
Other Names	ZNF563; zinc finger protein 563
Accession No.	Swiss-Prot#:Q8TA94NCBI Gene ID:147837
Calculated MW	Refer to figures
Formulation	Avoid freeze / thaw cycles.   Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

## Application Details

WB □ 1:500 - 1:2000

## Images



Western blot analysis of extracts of various cell lines, using ZNF563 antibody.

## Background

Note: This product is for in vitro research use only and is not intended for use in humans or animals.