

ZMIZ2 Polyclonal Antibody

Catalog No: #30011

Package Size: #30011-1 50ul #30011-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	ZMIZ2 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Hu,Ms,Rt
Immunogen Description	Recombinant fusion protein of human ZMIZ2 (NP_001287888).
Other Names	ZMIZ2; NET27; TRAFIP20; ZIMP7; hZIMP7; zinc finger MIZ-type containing 2
Accession No.	Swiss-Prot#:Q8NF64NCBI Gene ID:83637
Calculated MW	115kDa
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

Application Details

WB □ 1:500 - 1:2000

Images



Western blot analysis of extracts of various cell lines, using ZMIZ2 antibody.

Background

ZMIZ2 and ZMIZ1 (MIM 607159) are members of a PIAS (see MIM 603566)-like family of proteins that interact with nuclear hormone receptors. ZMIZ2 interacts with androgen receptor (AR; MIM 313700) and enhances AR-mediated transcription (Huang et al., 2005 [PubMed 16051670]).

Note: This product is for in vitro research use only and is not intended for use in humans or animals.