

## VPS37B Polyclonal Antibody

Catalog No: #30170

Package Size: #30170-1 50ul #30170-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)


## Description

Product Name	VPS37B Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Human,Mouse,Rat
Immunogen Description	Recombinant fusion protein of human VPS37B (NP_078943.1).
Accession No.	Swiss-Prot#:Q9H9H4NCBI Gene ID:79720
Calculated MW	31kDa
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

## Application Details

WB □ 1:500 - 1:2000

## Images



Western blot analysis of extracts of Rat thymus, using VPS37B antibody at 1:1000 dilution.

## Background

Note: This product is for in vitro research use only and is not intended for use in humans or animals.