

## DEPDC7 Polyclonal Antibody

Catalog No: #30187

Package Size: #30187-1 50ul #30187-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)


## Description

Product Name	DEPDC7 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Hu
Immunogen Description	Recombinant fusion protein of human DEPDC7 (NP_001070710.1).
Other Names	TR2, dJ85M6.4
Accession No.	Swiss-Prot#:Q96QD5NCBI Gene ID:91614
Calculated MW	58kDa
Formulation	Avoid freeze / thaw cycles.   Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

## Application Details

WB □ 1:500 - 1:2000

## Images



Western blot analysis of extracts of various cell lines, using DEPDC7 antibody.

## Background

Note: This product is for in vitro research use only and is not intended for use in humans or animals.