

TSPY3 Polyclonal Antibody

Catalog No: #30221

Package Size: #30221-1 50ul #30221-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	TSPY3 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Rat
Immunogen Description	Recombinant fusion protein of human TSPY3 (NP_001071165.2).
Other Names	CT78;TSPY3
Accession No.	Uniprot:P0CV98GeneID:728137
Calculated MW	35kDa
SDS-PAGE MW	35kDa
Formulation	PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.

Application Details

WB □ 1:500 - 1:2000

Images



Western blot analysis of extracts of Rat testis, using TSPY3 antibody.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.