

EGF Receptor (Phospho-Tyr998)

Catalog No: #12798



Package Size: #12798-1 50ul #12798-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	EGF Receptor (Phospho-Tyr998)
Purification	Affinity purification using immunogen.
Applications	WB
Species Reactivity	Hu Ms
Specificity	EGF Receptor (Phospho-Tyr998) detects endogenous levels of EGF Receptor only when phosphorylated at Tyr998
Immunogen Type	Peptide
Immunogen Description	Synthesized peptide derived from human EGF Receptor around the phosphorylation site of Tyr998.
Target Name	EGF Receptor (Phospho-Tyr998)
Modification	Phospho
Other Names	EGFR, EGFR receptor, ERBB, ERBB1, MENA, Panitumumab, PIG61, Proto-oncogene c-ErbB-1, V-erb-b homolog, EGF-R, HER1
Accession No.	Swiss-Prot#: P00533NCBI Gene ID: 1956
Target Species	human
Calculated MW	180 KD
Concentration	1.0mg/mL
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

Application Details

Western blotting: 1:1000-3000

Note: This product is for in vitro research use only and is not intended for use in humans or animals.